Micro infusion pump **EiPRE©IO®**

Implantable, Programmable and Refillable

Application Examples from peer reviewed publications



Elevate the best New Molecular Entity (NME) and New Chemical Entity (NCE) to the preclinical stage with the highest probability of successfully reaching the marketplace.

With a wide range of compatible solvents and vehicles, iPRECIO Micro Infusion Pumps are used as a platform in-vivo technology to:-


- Evaluate both non-optimized and optimized compounds.
- Maximize the information content about target engagement. (PK/PD)
 Minimize confounding variables due to stress and human interactions.
- Simulate humanized exposure profiles for translational research.



iPRECIO Micro Infusion Pumps

<Off-the-shelf> development tool for use in drug discovery

- 1. Reduced drug requirements
- 2. Large selection of compatible solvents used in drug discovery
- 3. Easy to modulate and time exposure profiles with non-optimized compound
- 4. Easy to use/program
- 5. Available since 2007



- Solubility issues and need a higher flow rate?
- Need more control for dosage due to narrow therapeutic index?
- Want to program a drug holiday or maximize efficacy/reduce toxicity with a timed dose during the mornings?
- Want to refill with a different test article/drug (sequential administration)?

1







The ability to program the device to start, stop and deliver different doses at different time points or just deliver one continuous dose makes iPRECIO ideally suited to the drug discovery and basic research process. All programmed in an easy to use PC based application software.

iPRECIO Micro Infusion Pumps for Drug Delivery

Exposure-enabling technology for advancing early preclinical studies and basic research

- Enables simple and complex dosing regimens at the click of the mouse/keyboard (ubiquitous PC) – several clicks
- · Automation which minimizes animal handling
- · Reduces stress and behavior anomalies
- · Parenteral route which is practical and extremely important

Basic requirements

- Surgical skills/training (important for successful use of iPRECIO Micro Infusion Pumps)
- Basic computer skills/literacy

Resources available from Primetech

- Surgical training videos and step by step Surgical Technical Notes
- User Manual, workflows and step by step programming guide
- Compatible vehicle/solvents and easy to use compatibility test kit

What researchers are saying:

Ease of programming: "I was pleasantly surprised with how easy it was to program, fill, and implant the pumps."

Programmable & implantable pump: "This device enables implementation of infusion protocols to reliably and precisely achieve the desired exposure profiles (shapes and timing) with low degree of invasiveness."

Improved drug delivery: "The infusion pumps enhanced the delivery of the drug and allowed for us to identify a clean behavioral antidepressant effect, devoid of complications due to daily injections."

Improved drug efficacy: "This study demonstrated that an equivalent effect was possible at a much lower dose than was previously studied (25µg/serotonin hydrochloride/kg/min) in the sham and DOCA–salt rat."

Reproducible results: "I have accumulated few more very nice recordings using iPrecio. Few recording are really breath-taking by reproducibility of responses."

"Your pump is AMAZING in terms of being able to do an intra-animal dose response curve. I absolutely, positively loved this. As a pharmacologist, there is nothing better."

Lead Optimization Study: "...we use them for studies to understand the PK-PD relationship of specific molecules. In terms of the infusion protocol it would be multiple steps to achieve a specific PK concentration in a PD study."

Things went well with the last iPRECIO study. The pumps did a fantastic job as they were programmed to do. iPRECIO data were in line with predicted/calculated values. As a matter of fact, we are in the process of completing another study using the iPRECIO pumps.

Contents

- 1. Introduction > P.1
- 2. What Researchers are saying > P.3
- 3. Contents > P.4
- 4. Research Applications 1-6 > P.5
 - a Liver-derived ketone bodies are necessary for food anticipation, Nature Communication
 - Dosing Profile Profoundly Influences Nicotinic Acid's Ability to Improve Metabolic Control in Rats,
 J. of Lipid Research
 - c The abruptness of terminating nicotinic acid delivery has a profound effect on free fatty acid and insulin rebound in rats, 51st EASD Annual Meeting, Stockholm 2015
 - d Ecto-domain phosphorylation promotes functional recovery from spinal cord injury, Scientific Reports
 - e Intrathecal administration using the iPRECIO® implanted pump
 - f Enhanced Resistance to Permeability Transition in Interfibrillar Cardiac Mitochondria in Dogs: Effects of Aging and Long Term Aldosterone Infusion, American Journal Of Physiology
 - g Highly Effective Auger-Electron Therapy in an Orthotopic Glioblastoma Xenegraft Model using Convection-Enhanced Delivery
- 5. GLP Studies with iPRECIO Pumps (Dual, SMP-200) > P.15
- 6. Toxicology Studies with iPRECIO Pumps (SMP-200) > P.15
- 7. Application Examples > P.16
- 8. Webinar : Gold Standard Physiological Measurements and Novel Drug Delivery Methods > P.17
- 9. The Ultimate Choice for Neuroscience > P.19
- 10. Recently in the news : > P.22
 - Miniaturized Neural System for Chronic, Local Intracerebral Drug Delivery (MiNDS)
 - <US Patent App.> 20180161282 Methods & compositions for the treatment of seizure-related disorders.
 - MitoTEMPO normalized ROS and prevented heart failure and arrhythmias
- 11. iPRECIO Micro Infusion Pumps for Cancer Research > P.26

Program what you need for overcoming:-

- a. Narrow therapeutic index
 - i. Maximizing efficacy with timed infusion
 - ii. Minimizing toxicity with timed infusions
 - iii. Flexibility to program 101 discreet infusion flow-rates. (SMP-300, 0.0 to 10.0µl/hours)
- 12. Selected Further Reading I: > P.27
 - a. Detailed Background iPRECIO and Drug Discovery
 - b. For More Trends
- 13. Selected Further Reading II: > P.28
 - c. Drug Discovery: iPRECIO Micro infusion pumps referenced
 - d. Drug Discovery PK/PD
- 14. Example Pump implantation site and drug administration site > P.29
- 15. Support Materials > P.29
 - a. Technical Note / Surgical Protocol
 - b. Surgical Videos
- 16. Product Information > P.31
- 17. Compatible solvents > backpage

Research Applications 1

SC Administration

Timed release of Test Article (TA), ßOHB

Rohit Chavan, Céline Feillet, Sara S. Fonseca Costa, James E. Delorme, Takashi Okabe, Jürgen A.Ripperger & Urs Albrecht

Liver-derived ketone bodies are necessary for food anticipation.

Nature Communications 7, Article number: 10580 doi:10.1038/ncomms10580

http://www.nature.com/ncomms/2016/160203/ncomms10580/full/ncomms10580.html?WT.ec_id=NCOMMS-20160205




Figure 1 (Figure 4a in Full Article) **Rescue of food anticipation in L Per2**-¹-**mice by ß-hydroxybutyrate**. (a) Timed release of ßOHB (green) but not NaCl (white) or Na-Pyruvate (purple) in L Per2-¹- mice mimics the ßOHB levels in plasma of L Per2+¹+ control animals (black). Measured after 15 days of infusion. Figure reproduced from Chavan et al. in Nature Communications as reference previously.

Figure 1 is reproduced from Liver-derived ketone bodies are necessary for food anticipation. http://www.nature.com/ncomms/2016/160203/ncomms10580/full/ncomms10580.html?WT.ec_id=NCOMMS-20160205 under Creative Common Attribution 4.0 International (CC BY 4.0) http://creativecommons.org/licenses/by/4.0/.

No changes were made for reproduction from Figure 4a of Chavan et al.

Purpose of the study:

Researchers were interested to know where Food Anticipation (FA) signals originate and what role components of the circadian clock might play. To test the potential of ßOHB as FA signal, iPRECIO SMP-300 programmable minimpumps were used to release ßOHB s.c. 6 hours prior to meal time under Restricted Feeding (RF) at ZT22 to reach a concentration normally observed in WT mice under RF preceeding feeding time.

iPRECIO SMP-300 pumps were used to test the potential of ßOHB as a FA signal.

Short methods or use of the pumps:

iPRECIO SMP-300 pumps were programmed to infuse saline vehicle at 2 ul/h, or D-\(\text{SOHB} \) at 2 ul/h, or Sodium pyruvate at 5 ul/h, or coconut oil at 5 ul/h prior to meal time (6 h, ZT22-ZT4) under Restricted Feeding (RF)

Results/significance:

Liver-derived ketone bodies are necessary for food anticipation.

Timed Release of &OHB partially rescues FA.

Research Need:

Timed Release of ßOHB in free moving animal with minimum or no handling to reduce stress and any confounding effects.

Additional information on mini-pump implant

Male and female L Per2^{+/+} and L Per2^{-/-} mice (3-5 months old)Telemetry transmitter (G2 Emitter) was i.p. implanted in each mouse under gaseous anaesthesia. At least 10 days after the transmitter implantation an iPRECIO programmable micro infusion pump (SMP/UCD 300; Primetech Corp., Japan) was implanted in subgluteal space(s.c. administration) on the back of each L Per2^{-/-} mouse. Subcutaneous administration.

Related Circadian rhythm Research using iPRECIO SMP-200 in mice

In vivo imaging of clock gene expression in multiple tissues of freely moving mice Nature Communications 7, Article number: 11705 doi:10.1038/ncomms11705

https://www.nature.com/articles/ncomms11705

Research Applications 2

SC Administration

Can different dosing (12 hour rectangular exposure profile) and terminating profile (a slow-step down) of Nicotinic Acid (NiAc) prevent/delay tolerance development and attenuate the FFA rebound development respectively.

Tobias Kroon (2016) PhD Thesis,

<Optimizing Nicotinic Acid Delivery for Durable Anti-lipolysis and Improved Metabolic Control>,

http://pub.epsilon.slu.se/13324/

http://pub.epsilon.slu.se/13324/1/kroon_t_160429.pdf

Thesis and publications cover Drug Discovery implications

- 1. Importance of time-series disease model
- 2. Continuous vs. intermittent drug exposures /Programmable, implantable mini-pump
- 3. Time exposure to physiology/Shape of exposure
- 4. Meta-analysis/Rank candidates/Predict designs

Tobias Kroon, Ann Kjellstedt, Pia Thalén, Johan Gabrielsson, Nicholas D. Oakes

Dosing Profile Profoundly Influences Nicotinic Acid's Ability to Improve Metabolic Control in Rats

The Journal of Lipid Research, doi: 10.1194/jlr.M058149 , July 13, 2015

http://www.jlr.org/content/early/2015/07/13/jlr.M058149.abstract

Kroon T, Baccega T2, Olsén A, Gabrielsson J, Oakes ND

Nicotinic acid timed to feeding reverses tissue lipid accumulation and improves glucose control in obese Zucker rats [S].

J Lipid Res. 2017 Jan; 58 (1): 31-41 Doi: 10.1194 / jlr.M 068395. Epub 2016 Nov 15.

https://www.ncbi.nlm.nih.gov/pubmed/27875257






Fig. 2 (Figure 1 in Kroon et al.)
A: NiAc and saline infusion profiles across studies I–III.

Black (NiAc) and open (saline) bars represent time periods of constant rate infusions during days 1–5. B: Terminal protocol for studies I (NiAc-induced FFA lowering) and III (NiAc-induced changes in adipose tissue gene expression).

C: Terminal protocol for study II (hyperinsulinemic-isoglycemicla clamps).




Fig. 3. (Figure 5 in Kroon et al.) Plasma insulin (A, B) and glucose (C. D) concentration in lean (left) and obese (right) following infusion of saline (lean n = 5, obese n = 12) or NiAc (0.17µmol·min -1·kg -1) given acutely (NiAc naïve, n = 7/group) or following 5 days continuous (Cont. NiAc, lean n = 4, obese n = 8) or intermittent (Inter. NiAc, lean n = 4, obese n = 9) or 11 days intermittent (Inter. NiAc Day 11, obese n = 4) dosing. The black horizontal bar represents the period of acute NiAc/saline infusion. Data presented as mean ± SE.

Figures 2 and 3 licensed material. © <2015> The American Society for Biochemistry and Molecular Biology. Warranties: None Publisher makes no representations or warranties with respect to the licensed material and adopt on its own behalf the limitations and disclaimers established by CCC on its behalf in its Billing and Payment terms and conditions for this licensing transaction.

Purpose of the study:

Dosing Profile Profoundly Influences Nicotinic Acid's Ability to Improve Metabolic Control in Rats

Researchers wanted to compare the ability of continuous versus intermittent NiAc administration to suppress FFA levels in metabolic healthy and insulin-resistant rats.

The abruptness of terminating nicotinic acid delivery has a profound effect on free fatty acid and insulin rebound in rats

The aim of this study was to determine whether a slow step-down NiAc infusion protocol (Step-Down group) vs. simply turning infusion off (On/Off group) could attenuate the FFA rebound development.

iPRECIO SMP-200 pumps were programmed to deliver the required exposure profiles of Nicotinic Acid to study impact on tolerance development (see figure 2A) and attenuate the FFA rebound development respectively (not shown)

Results/significance:

An Intermittent NicAc dosing strategy succeeded in retaining FFA lowering and improving insulin sensitivity in obese Zucker rats. Gradual step-down reduction of NiAc infusion actually degraded the anti-lipolytic effectiveness of NiAc compared to abrupt withdrawal.

Research Need:

Ability to quickly and easily adjust dosing profiles based on PK and PD effects and deliver doses without stressors which could change metabolic activity of animals.

Intrathecal Administration

Research Applications 3

Continuous infusion of PKA and ATP at 1µl/hour for 14 days where solution in pump was changed every 2 days due to stability of PKA and ATP.

Kenji Suehiro, Yuka Nakamura, Shuai Xu, Youichi Uda, Takafumi Matsumura, Yoshiaki Yamaguchi, Hitoshi Okamura, Toshihide Yamashita & Yoshinori Takei

Ecto-domain phosphorylation promotes functional recovery from spinal cord injury

Scientific Reports 4, Article number: 4972 (2014) doi:10.1038/srep04972

http://www.nature.com/articles/srep04972




Figure 4 is reproduced from Ecto-domain phosphorylation promotes functional recovery from spinal cord injury

http://www.nature.com/articles/srep04972

under Creative Common Attribution 4.0 International (CC BY 4.0) http://creativecommons.org/licenses/by/4.0/.

No changes were made for reproduction from Fig. 1 of Suehiro et al.

Figure 4 (Figure 1 from Suehiro et al)

I Treatment with PKA plus ATP diminishes damage from traumatic SCI. (a) The depth of injury and location of sections used in (c) are illustrated schematically. The dorsal corticospinal tract (dCST) and the dorsolateral corticospinal tract (dICST) were severed. (b) The BBB scores of vehicle-treated, PKAtreated, ATP-treated and PKA1ATPtreated SCI rats were assessed at the indicated days after SCI. The points on the graph indicate the average BBB score from six independent rats, and the error bars indicate the standard deviation (S.D.) (*p, 0.05, **p, 0.01, ***p, 0.001 vs.vehicle-treated rats, Student's t-test). (c) The BDAlabelled dCST was visualised. Images are taken from transverse sections at either 5 mm caudal or rostral to the lesion, as shown in (a). The bar indicates 25 mm. (d) The number of BDA-positive axons at T8 or T10 was normalised to the number of BDA positive axons at C1 (intact region of the spinal cord). The average and the S.D. from three independent animals are shown. No significant differences between the vehicle-treated rats and the PKA/ATP-treated rats were observed (*p. 0.05, **p. 0.01,

Purpose of the study:

Investigate if inhibition of Nogo-66 receptor (NgR) via ecto-domain phosphorylation by protein kinase A (PKA), which blocks activation of the receptor can promote recovery following spinal cord injury.

iPRECIO SMP-200 pumps were used to infuse PKA plus ATP for 14 days at 1µl/hour. Solution in reservoir was changed every 2 days.

Results/significance:

Authors found that infusion of PKA plus ATP into the damaged spinal cord can promote recovery of locomotor function.

Research Need:

Ability to replace unstable test articles or drugs easily and rapidly without additional surgeries and stress.

Related publication examples: Refilling to Improve Test Article Stability

Hemoglobin induced lung vascular oxidation, inflammation, and remodeling contributes to the progression of hypoxic pulmonary hypertension and is attenuated in rats with repeat dose haptoglobin administration

Free Radical Biology and Medicine D Irwin et

al.doi:10.1016/j.freeradbiomed.2015.01.012

http://www.sciencedirect.com/science/article/pii/S0891584915000192

Free hemoglobin induction of pulmonary vascular disease: evidence for an inflammatory mechanism.

Am J Physiol Lung Cell Mol Physiol. 2012 Aug;303(4):L312-26. Epub 2012 Jun 22.

http://www.ncbi.nlm.nih.gov/pubmed/22728465

Research Applications 4

Intrathecal Administration

Excerpt from Mitchell et al. Full reference in box.

Regulatory request to perform an epidural and/or intrathecal animal study to assess degradents associated with a pharmaceutical product that was given epidurally in humans.

Mitchell D., Read, K., Chapman M. and Patten D.

Intrathecal administration using the iPRECIO® implanted pump

Development in Life Sciences, Vol 14, No. 4

http://cdn2.hubspot.net/hubfs/212573/docs/Envigo/Envigo Pharma Dils 14.4.4.pdf?t=1460116975327

Purpose of the study:

The customer requested a rat study involving intrathecal infusion for 72-hours of two different degradent mixtures and appropriate controls with acute and delayed endpoints and investigations of local and systemic toxicity. Clinical relevant concentrations of degradents to attain comparable exposure with humans would be necessary.

Developments in LifeSciences Vol. 14 No. 4

Intrathecal administration using the iPRECIO[®] implanted pump

David Mitchell BSc (Hons) DABT, Senior Toxicologist, Toxicology Operations, Envigo, UK.
Kate Read MA VetMB MRCVS, Veterinary Clinician, Veterinary Services, Envigo, UK.
Melissa Chapman BSc (Hons), Senior Study Director, Toxicology Operations, Envigo, UK.
Duncan Patten FIAT RAnTech, Associate Director, Laboratory Animal Technologies, Envigo, UK.

The background to this project was a regulatory request to perform an epidural and/or intrathecal animal study to assess degradents associated with a pharmaceutical product that was given epidurally in humans. There was a concern that there might be inadvertent intrathecal administration of the product and degradents. The customer requested a rat study involving intrathecal infusion for 72-hours of two different degradent mixtures and appropriate controls, with acute and delayed toxicity endpoints and investigations of local and systemic toxicity. We had

could compromise welfare, in particular for clinical signs associated with increased intrathecal pressure.

The optimal solution was to use the

iPRECIO* SMP-200 programmable

peristaltic pump implanted

subcutaneously and linked to an

intrathecal catheter

iPRECIO SMP-200 pumps were used to infuse 1µl/hr of artificial CSF intrathecally following surgery and during the recovery period. Animals recovered well with no adverse clinical signs in the post –operative period. During the treatment period; infusion at 30µl/hr, a small number of animals (5 out of 72) showed hindlimb paresis. Examination of aspirated dose volumes demonstrated accurate pump function.

Results/significance

This method (iPRECIO SMP-200 linked to an intrathecal catheter) is suitable for controlled continuous infusion into the intrathecal space of the rat. The surgical procedure is reproducible and considered to be less invasive than intrathecal access via the cisterna magna. The use of the programmable iPRECIO® pump allows for an ambulatory infusion model without the need to tether the animals. This permits behavioural assessment and is an improvement in animal welfare; animals are able to display normal behaviours post operatively.

Research Need:

A standard method for intrathecal infusion in industry and academia which would not be a confounding factor in the assessment of CNS endpoints (modified Irwin assessment). The infusion system must provide a suitable flow rate over at least 72 hours.

- •The pump must allow the flexibility to start infusion immediately following surgery or at a later time.
- •The pump must have a reservoir that can be evacuated and refilled, percutaneously, by syringe and needle so there would be the opportunity for a period of recovery from surgery before administration of the degradant mixtures while avoiding the risk of catheter occlusion by administering saline or artificial cerebrospinal fluid.

Research Applications 5

Jugular Vein (IV) Administration

Aldosterone was continuously infused with SMP-200 programmable infusion pump that delivered aldosterone into the jugular vein. D-Aldosterone was infused into the jugular vein at a dose of 30 μ g·kg⁻¹·day⁻¹ in a solution of 15% ethanol,

50% DMSO, and 35% water at a concentration of 10 mg aldosterone/ml.

Enhanced Resistance to Permeability Transition in Interfibrillar Cardiac Mitochondria in Dogs: Effects of Aging and Long Term Aldosterone Infusion.

Am J Physiol Heart Circ Physiol aipheart.00674.2012;

http://ajpheart.physiology.org/content/early/2012/12/10/ajpheart.00674.2012.abstract?sid=1c8187a4-b1a5-41e2-9e88-20610af15128

Purpose of the study:

Effect of aging and long-term aldosterone infusion on respiratory function and resistance to mitochondrial permeability transition (MPT) in subsarcolemmal and interfibrillar cardiac mitochondria (SSM and IFM) from healthy young (1 year) and old (8 year) female beagles.

iPRECIO SMP-200 pumps were used to infuse Aldosterone for 14 weeks at a dose of 30 μg·kg⁻¹·day⁻¹ The pump reservoir was 900 μl and was refilled percutaneously every 20–30 days through an injection port on the pump. The pump reservoir was evacuated before refilling to ensure the pump had properly discharged its contents and was then refilled using a 26-gauge needle. This procedure was done in conscious animals with no evidence of discomfort.

Results/significance

Authors demonstrated in a large animal model that resistance to MPT is greater in IFM than in SSM in young and old female dogs. When old dogs were stressed with aldosterone infusion, there was selective enlargement of SSM and greater susceptibility to MPT, with no change to IFM.

Research Need:

Long term/chronic 14 week infusions with the ability to refill and check performance of implanted pumps.

Research Applications 6

Brain Administration

Pumps were programmed to instant mode, constant mode and 5µl/hour infusion rate.

They were initially loaded with isotone saline or 0.1 mM MTX. Two days later, residual saline or MTX was extracted from the pump reservoirs and refilled with 960 μ l of 0.3 μ g/ml ¹²⁵l-UdR or ¹²⁷l-UdR. See figure 6 below for results obtained. Reproduced with permission from Thisgaard et al. (CC BY-NC-ND 4.0).

Thisgaard et al.

Highly Effective Auger-Electron Therapy in an Orthotopic Glioblastoma Xenograft Model using Convection-Enhanced DeliverHighly Effective Auger-Electron Therapy in an Orthotopic Glioblastoma Xenograft Model using Convection-Enhanced Deliver

Theranostics 2016, Vol. 6, Issue 12 2016; 6(12): 2278-2291. doi: 10.7150/thno.15898

http://www.thno.org/v06p2278.htm




Figure 6. Kaplan-Meier plot showing that the survival benefit of neoadjuvant MTX + 125|-UdR as stand-alone Auger-therapy (group4) or with concomitant, systemic chemotherapy (group5) was highly significant compared with the non-radioactive, but chemically identical treatment MTX + 127I-UdR (group3, p=0.0001 and p<0.0001, respectively) or untreated controls (group1, both p<0.0001). The Auger-therapy was also significantly better than systemic TMZchemotherapy alone (group6, p=0.0001). Reproduced with permission from Thisgaard et al. (CC BY-NC-ND 4.0).

Purpose of the study:

The overall aim of this was to test the effect and safety profile of ¹²⁵I-UdR therapy in vitro and in vivo on immature Glioblastomas (GBMs) spheroid cultures (GSCs) and orthotopic xenografted GBM-bearing rats, respectively. A further objective was to determine if further therapeutic effect was achieved when combining ¹²⁵I-UdR therapy with the currently used first-line chemotherapeutic agent TMZ.

Pumps were initially loaded with isotone saline or 0.1 mM MTX. Two days later, residual saline or MTX was extracted from the pump reservoirs and refilled with 960 μ l of 0.3 μ g/ml ¹²⁵I-UdR or ¹²⁷I-UdR.

Results/significance

The multidrug approach including CED of MTX and the AEE-compound ¹²⁵I-UdR in combination with systematic TMZ was safe and very effective in the orthotopic xenograft GBM model, leading to 100% survival.

Research Need:

The ability to evaluate combinational therapy/multidrug approach easily and rapidly without additional surgeries and stress.

GLP Studies with iPRECIO Pumps

Laura Ringer

The use of the iPRECIO Dual Inlet Infusion Pump in Ambulatory Cardiovascular Dog Studies

DSI East Coast User Group Meeting, Philadelphia, PA, United States October 29th and 30th 2015 The Use of the Iprecio Dual Inlet Infusion Pump in Ambulatory Cardiovascular Dog Studies

Laura Ringer
Pizer

Yearn Members. Pricer Harrs, Virgent Bernardo


Duncan Patten (Huntingdon Life Sciences, UK)

Use of iPRECIO implantable micro infusion pumps in rats

4th Infusion Technology Organization Meeting, May 8th-9th 2014, Harrogate, UK.

Perron J., Frenette V., and Copeman C.

Validation and use of the iPRECIO[®] Micro Infusion Pump on GLP studies Society of Toxicology Annual Meeting, San Francisco, United States, March 11th to 14th 2012.



Toxicology Studies with iPRECIO Pumps

Masaru Tsuboi, Yoshihide Ueda, Yasufumi Ota, Hiroshi Takehara, Takuya Aoshima, Fukutaro Mizuhashi


Physiological conditions in iPRECIO® -implanted rats

Fundamental Toxicological Sciences Vol.3 (2016) No.1 p.1-8

https://www.jstage.jst.go.jp/article/fts/3/1/3_1/_article

Application Examples

5-HT dose response with control period : 5-25 greater sensitivity



Drug Delivery: Enabling Technology for Drug Discovery and Development.

iPRECIO® Micro Infusion Pump:


Programmable, Refillable, and Implantable

Tsung Tan, Stephanie W. Watts, and Robert Patrick Davis Front Pharmacol. 2011; 2: 44. Published online 2011 July 29.

doi: 10.3389/fphar. 2011.00044

https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3149148/

Dose response: Dobutamine, verapamil & saline 3 test articles per animal (pump)




Drug Delivery: Enabling Technology for Drug Discovery and Development. iPRECIO® Micro Infusion Pump:

Programmable, Refillable, and Implantable

Tsung Tan, Stephanie W. Watts, and Robert Patrick Davis Front Pharmacol. 2011; 2: 44. Published online 2011 July 29. doi: 10.3389/fphar. 2011.00044

https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3149148/

100nl bicuculline methiodide (BMI) bolus injections



Zaretsky D.V., Zaretskaia M.V., Durant P.J., Rusyniak D.E.
The use of microinfusion pump to perform intrahypothalamic injections in conscious

Neuroscience 2012, New Orleans, USA., October 13th - 17 2012


http://www.abstractsonline.com/Plan/ViewAbstract.aspx?mID=2964&sKey=87d8b951-316f-466a-9eb7-4b154d0bbd2c&cKey=b4b8338f-9bd2-44e2-bcf4-6e05a36cbbcb&mKey=%7b70007181-01C9-4DE9-A0A2-EEBFA14CD9F1%7d

Comparison of arterial pressure and plasma ANG II responses to three methods of subcutaneous ANG II administration

Comparison of arterial pressure and plasma AngII responses to three methods of subcutaneous AngII administration Kuroki M.T. . , Gregory D. Fink , John W. Osborn

American Journal of Physiology - Heart and Circulatory PhysiologyJul 2014, DOI: 10.1152/ajpheart.00922.2013 http://ajpheart.physiology.org/content/early/2014/06/30/ajpheart.00922.2013

Webinar













the Q&A Report

the iPRECIO Resources

Learn how iPRECIO Micro Infusion Pumps will help you with your research needs. Webinar recordings, slide decks, iPRECIO resources (Bibliography of Scientific publications by product features) and Q&A report are now all < Open Access >





Dr Christian Schnell

Want to see objective data which lead Dr Schnell to the following conclusions?

- 1. The best situation is to get rid of stress because training is not really a solution.
- An implantable pump is the only way to deliver the compound without interfering at the moment of delivery.

The open access on-demand webinar, resources, and Q&A are available 24/7.



Dr Robert Doyle

How are we going to know that if we do SC or Oral administration that we are not inducing a stress response which essentially negates the data we collect at the back end in terms of food intake and body weight as they pertain to our study?

Q: What were the best resources for learning how to program the iPRECIO mini pumps? Was it relatively simple?

Ans. [R. Doyle]:

Highly skilled technician used iPRECIO for the first time and did not have any difficulties and expect that similar skilled technicians to be fine. Medicinal chemists with only basic husbandry have been trained by technician (Monday–Friday week), to program and implant (surgery). Never lost a pump and never had a student screw-up.

What drug delivery concerns did Dr. Robert Doyle have and how did he solve it? The open access on-demand webinar, slide deck, iPRECIO resources and Q&A are available 24/7. See also publication in Endocrinology.

- 1. We have to be particularly conscious of the fact that when looking at <Food Intake> and <Body Weight Reduction>, the handling of the animals is particularly confounding because inducing stress or inducing fight-or-flight response etc. is going to play havoc with our data collection.
- 2. And if we are dealing with subtle changes that can be the reason we lose out on seeing something we wouldn't otherwise and obviously that's something we were really really conscious of and concerned about when we began to do this project.
- 3. How are we going to know that if we do SC or Oral administration that we are not inducing a stress response which essentially negates the data we collect at the back end in terms of food intake and body weight as they pertain to our study?

Endocrinology Publication

Vitamin B₁₂ Conjugation of Peptide-YY₃₋₃₆ Decreases Food Intake Compared to Native Peptide-YY₃₋₃₆ Upon Subcutaneous Administration in Male Rats

Endocrinology February 6, 2015 doi: 10.1210/en.2014-1825

https://doi.org/10.1210/en.2014-1825 [Open Access]

SMP-200 iPRECIO Use

- ·Recovery period after surgery of 7 days and 5 day baseline period (2µl/hr of saline)
- Treatments were delivered sc with five pulses per day; three 1 hour pulses of 10 nmol.kg⁻¹.h⁻¹ (20 μl.h⁻¹) with three hours between pulses and two 1 hour pulses of 5 nmol.kg⁻¹.h⁻¹ (10 μl.h⁻¹) with five hours between pulses.

The Ultimate Choice for Neuroscience

iPRECIO Micro Infusion Pumps for Drug Delivery Implantable Programmable Refillable

- ·The only way to deliver compound without interfering at the moment of delivery
- · Paired data sets: Program a recovery/baseline period prior to drug delivery for control period for comparison.
 - > Recovery period after surgery (pump stop or saline infusion)
 - > Baseline period (pump stop or saline infusion)
 - > Drug delivery /Treatment period (start pump or exchange from saline to drug)
 - Continuous
 - Intermittent
 - Dose escalation / de-escalation
 - Circadian
 - > Reversibility (pump stop or exchange to saline)
- · Infuse directly to brain
- · Infuse directly to intrathecal space
- ·SC, IP and IV administration

Example Drug Delivery Regimen (Figure 1 reproduced from Thisgaard et al. (CC BY-NC-ND 4.0) Schedule what you require: program and/or exchange infusate as per study requirements




Figure 1 reproduced from Thisgaard et al. (CC BY-NC-ND 4.0)

Highly Effective Auger-Electron Therapy in an Orthotopic Glioblastoma Xenograft Model using Convection-Enhanced Delivery Thisgaard et al. Theranostics 2016, Vol. 6, Issue 12 2016; 6(12): 2278-2291. doi: 10.7150/thno.15898

http://www.thno.org/v06p2278.pdf Attribution-Non Commercial-No Derivatives 4.0 International (CC BY-NC-ND 4.0) https://creativecommons.org/licenses/by-nc-nd/4.0/

Selected CNS Publications

Intrastriatal Memantine Infusion Dampens Levodopa-Induced Dyskinesia and Motor Deficits in a Mouse Model of Hemiparkinsonism

BRIEF RESEARCH REPORT ARTICLE Front. Neurol., 05 December 2019

https://www.frontiersin.org/articles/10.3389/fneur.2019.01258/full

Key words: intracerebral brain infusion, levodopa-induced dyskinesia, memantine, N-methyl-D-aspartate receptor, Parkinson's disease





Figure 1A reproduced without modification from doi.: 10.3389/fineur.2019.01258. Attribution 4.0 International (CC BY 4.0) https://creativecommons.org/licenses/by/4.0/

Continuous cerebroventricular administration of dopamine: A new treatment for severe dyskinesia in Parkinson's disease?

Neurobiology of Disease, Vol. 103, 2017, 24-31

http://dx.doi.org/10.1016/j.nbd.2017.03.013

Pump setting delivery in 6-OHDA rats:



Supplementary Image 1.

Dopamine delivery from the pump trough the rat brain cannula begin each day from zeitgeber time -10h (i.e. 9pm) to zeitgeber time 6 h (i.e. 1pm), over 16h during 30 days.

Image 1 & text reproduced without modification from C. Laloux et al. (CC BY-NC-ND 4.0) https://creativecommons.org/licenses/by-nc-nd/4.0/

Rhythmic Release of Corticosterone Induces Circadian Clock Gene Expression in the Cerebellum

Neuroendocrinology. 2019 Sep 27. doi: 10.1159/000503720

https://www.karger.com/Article/Abstract/503720

Key words: Cerebellum, Circadian, Clock gene, Corticosterone, Suprachiasmatic nucleus,iPRECIO programmable micropump

Investigate drug-evoked adaptations with different patterns of exposure.

Selected CNS Applications with iPRECIO Micro Infusion Pumps

Addiction/ Drug abuse liability

- · Adversive effects of drug withdrawal in rats and mice
- · Withdrawal Test
 - > Test potential compounds which may have similar effects in the same animals or reduce the signs of withdrawal
 - > Abrupt cessation

Interruption of continuous opioid exposure exacerbates drug-evoked adaptations in the mesolimbic dopamine system

Neuropsychopharmacology (2020) | Published: 20 February 2020

https://doi.org/10.1038/s41386-020-0643-x Subjects: Addiction, Reward

Discrimination Learning in Oxycodone-Treated Nonhuman Primates

Drug and Alcohol Dependence, Available online 27 November 2019, 107778

https://doi.org/10.1016/j.drugalcdep.2019.107778

Keywords: Opioid, Oxycodone, Naltrexone, Self-administration, Withdrawal, Cognition, Nonhuman primate

Convergent and Divergent Behavioral Changes Caused by Different Patterns of Morphine Exposure in Mice

International Narcotics Research Conference (INRC), Chicago, 9 - 14 of July 2017 https://www.fourway.es/view/324/abstracts/#5203

Mechanisms specific to methamphetamine-associated memory disruption by nonmuscle myosin II inhibition

420.01 / SS23, Neuroscience Annual Meeting 2017, Nov 11th to 15th 2017 Washington DC http://www.abstractsonline.com/pp8/#J/4376/presentation/16699

CDKL5 PROTEIN SUBSTITUTION THERAPY RESCUES NEUROLOGICAL PHENOTYPES OF A MOUSE MODEL OF CDKL5 DISORDER

Human Molecular Genetics, ddy064, https://doi.org/10.1093/hmg/ddy064

https://academic.oup.com/hmg/advance-article-abstract/doi/10.1093/hmg/ddy064/4892297?redirectedFrom=fulltext

Recently in the news: MiNDS

Fortunately for us, iPRECIO® too. Playing our small part for Science.



CC BY-NC-ND*,
Credit M. Scott Brauer
Miniaturized Neural System for Chronic,
Local Intracerebral Drug Delivery
(MiNDS)

*CC BY-NC-ND: Attribution-NonCommercial-NoDerivatives 4.0 International

https://creativecommons.org/licenses/by-nc-nd/4.0/legalcode

https://www.media.mit.edu/projects/miniaturized-neural-system-for-chronic-local-intracerebral-drug-delivery/press-kit/

Miniaturized neural system for chronic, local intracerebral drug delivery

Science Translational Medicine 24 Jan 2018: Vol. 10, Issue 425, eaan2742

DOI: 10.1126/scitransImed.aan2742

http://stm.sciencemag.org/content/10/425/eaan2742

Focal, remote-controlled, chronic chemical modulation of brain microstructures

PNAS July 10, 2018 115 (28) 7254-7259; https://doi.org/10.1073/pnas.1804372115

Selected Patents: iPRECIO® Micro Infusion Pumps have been used/referenced.

Document	Document Title
US20180264191	SYSTEMS AND METHODS FOR NEURAL DRUG DELIVERY AND MODULATION OF BRAIN ACTIVITY
	A neural drug delivery system is disclosed. In an embodiment, the system includes two or more microtubes, each having a distal end, a proximal end, and elongate channel body extending http://www.freepatentsonline.com/y2018/0264191.html
US20170151193	PHARMACEUTICAL SOLUTION COMPRISING DOPAMINE FOR USE IN TREATING PARKINSON'S DISEASE The present invention is directed to pharmaceutical solution comprising at least dopamine for use in treating Parkinson's disease, wherein said pharmaceutical solution is kept under anaerobic http://www.freepatentsonline.com/y2017/0151193.html
US20170246200	MICRORNA-132/212 FOR THE TREATMENT OF NEURODEGENERATIVE DISORDERS The present relates to the use of miRNA-132/212 mimics or activators comprising a doubled stranded ribonucleic acid molecule comprising a seed region sequence of miRNA-132 or miRNA-212, a spacer, http://www.freepatentsonline.com/y2017/0246200.html

Recently in the news: United States Patent Application 20180161282

Methods and compositions for the treatment of seizure-related disorders by Adamas Biopharmaceutical Inc. (https://www.adamaspharma.com/) includes references to iPRECIO Micro Infusion Pump use in Examples 13 to 16. In example 13, "The PK parameters (Cmax, Tmax, T1/2, AUC) of 200 mg IR lacosamide BID and 400 mg ER lacosamide QD in humans can be replicated in rats using programmable subcutaneous infusion pumps.

- CNS side effects were compared when delivered as a continuous infusion in rats to simulate an extended release (ER) plasma profile or as a bolus to simulate an immediate release (IR) plasma profile.
- CNS side effects were measured as an impairment of motor coordination in rodents using the rotor-rod.

Excerpts from US Patent Publication 20180161282

Example 13: Pharmacokinetic Study of Lacosamide Administered by Continuous Subcutaneous Infusion to Rats for 7 Days

The objective of this study is to demonstrate that the PK profile from a standard regimen of 200 mg IR lacosamide BID in humans, along with a novel PK profile of 400 mg ER lacosamide QD, can be replicated in rats using programmable subcutaneous infusion pumps. Jugular vein cannulated male and female Sprague Dawley rats weighing between 240-260 grams are used for this study. Rats are implanted with programmable infusion pumps (iPRECIO SMP-200)

Example 14: - related to training for Rotor-rod Training and Test

Following the period of acclimation, rats are surgically implanted with subcutaneous programmable pumps (iPRECIO), per the manufacturer's instructions. Prior to implantation, pumps are pre-programmed and filled with sterile 0.9% saline according to the manufacturer's instructions. Implantation is done under anesthesia by using aseptic techniques.

Three days following pump implantation surgery, rats are pre-trained on the rotarod and assigned to treatment groups based on their body weights.

Infusion of lacosamide commences on day four post-pump implantation surgery. Baseline measurements are made on the rotarod on Day 4. Immediately thereafter, saline in the iPRECIO pumps is withdrawn and the pump is refilled with lacosamide solution (25 mg/mL in 0.9% saline). Pumps are refilled during the 7 day period as necessary, through the access port of the iPRECIO pump.

Excerpts from US Patent Publication 20180161282 continued.

for Simulated BID dosing

TABLE EX 13A

Study Design and Infusion Protocol

Group	Number of animals	Dose Period ²	Dose rate/rat ¹ (mg/h/rat)	Concentration in pump (mg/mL)	Pump Flow rate (µL/h)	Plasma Collection Times (hours)
1 4	0 - 2 h	0.59	20	29.5	1, 4, 8, 14, 20	
	4	2 - 12 h	0.12	20	6	84, 144, 158
2	4	0 - 2 h	0.59	20	29.5	2, 6, 12, 18
		2 - 12 h	0.12	20	6	24, 98, 156, 168

¹ Lacosamide dose estimates are based on average animal body weight of 250 g

Simulated ER QD Dosing

TABLE EX 13B

Study Design and Infusion Protocol

Number		Dose	Dose rate/rat ¹	Concentration	Pump Flow rate	Plasma Collection
Group of animals	Period ²	(mg/h/rat)	in pump (mg/mL)	(µL/h)	Times (hours)	
1	4	0 - 12 h	0.30	20	15	3, 9, 15, 21, 30
		4	12 - 24 h	0.15	20	7.5
2	4	0 - 12 h	0.30	20	15	6, 12, 18, 24, 36
		12 - 24 h	0.15	20	7.5	48, 144, 168

¹Lacosamide dose estimates are based on average animal body weight of 250 g

For more information, see https://patents.google.com/ or http://appft.uspto.gov.

² Pump infusion protocol is repeated every 12 hours for a total of 7 days


² Pump infusion protocol is repeated every 12 hours for a total of 7 days

Recently in the news: MitoTEMPO normalized ROS and prevented heart failure and arrhythmias in guinea pig model.

Fortunately for us, iPRECIO® too. Playing our small part for Science.

Enabling discovery and delivery.

Heart rate recordings over a period of 96 hours: showing internal consistency of iPrecio pumps Pumps inject isoproterenol once daily for a period of 1 hour. Heart rate recovers to baseline in 3-4 hours.



Full data and research published in

Mitochondrial ROS Drive Sudden Cardiac Death and Chronic Proteome Remodeling in Heart Failure

Circulation Research. 2018; CIRCRESAHA.118.312708, https://doi.org/10.1161/CIRCRESAHA.118.312708

"The iPRECIO pump enables reliable, timed delivery of pharmacological agents. The pump was programmed for 1 hour delivery of Isoproterenol at 30ul/hr, total dose of 2mg/kg/day, once a day at the same time (1pm)"

Related Reference.

Intracerebroventricular infusion of donepezil prevent cardiac remodeling and improves the prognosis of chronic heart failure rats

The Journal of Physiological Sciences volume 70, Article number: 11 (2020)

https://doi.org/10.1186/s12576-020-00739-0

iPRECIO Micro Infusion Pumps for Cancer Research

Solubility and Precipitation issues?

More difficult to dose correctly and need more control?

Program what you require

- · Solubility issues and need a higher infusion flow-rate to reduce drug concentration and precipitation risk
- · Difficult to dose correctly and need to be able to have accurate flow-rates/dose groups
- · Suited for intermittent dosing of onco substances – daily for 1 hour or every 2 days for 2 hours.
- ·Would like to allow tumor size to grow to a certain size before drug infusion
- · Want to program a drug holiday
- · Want to evaluate chrono release for maximum efficacy and minimize toxicity

Cancer Research Publications

Establishment of an orthotopic bladder cancer model to evaluate continuous intravesical delivery of small molecule inhibitors in the nude rat

AACR 106th Annual Meeting 2015; April 18-22, 2015; Philadelphia, PA

http://cancerres.aacrjournals.org/content/75/15_Supplement/5146.short

Convection-enhanced delivery of an anti-miR is well-tolerated, preserves anti-miR stability and causes efficient target de-repression: a proof of concept.

Journal of Neuro-Oncology 2015 Oct 1.

http://link.springer.com/article/10.1007%2Fs11060-015-1947-2

http://www.ncbi.nlm.nih.gov/pubmed/26428358



Tajiri et al. (Kyushu University, Japan)

Targeting Ras-Driven Cancer Cell Survival and Invasion through Selective Inhibition of DOCK 1

Cell Reports 19, 969-980, May 2, 2017

http://dx.doi.org/10.1016/j.celrep.2017.04.016

Reproduced from Tajiri et al. (CC BY-NC-ND 4.0) without modification https://creativecommons.org/licenses/by-nc-nd/4.0/ Conscious Talian Protein Racal Racal

Maxim Shevtsov et al.

Granzyme B Functionalized Nanoparticles Targeting Membrane Hsp70 - Positive Tumors for Multimodal Cancer Theranostics

Small, 2019 - Wiley Online Library

https://onlinelibrary.wiley.com/doi/abs/10.1002/smll.201900205

Selected Further Reading I

> For More detailed Background iPRECIO and Drug Discovery, see

Drug Delivery: Enabling Technology for Drug Discovery and Development. iPRECIO[®] Micro Infusion Pump: Programmable, Refillable, and Implantable

Front Pharmacol. 2011; 2: 44. Published online 2011 July 29. doi:10.3389/fphar.2011.00044

http://www.ncbi.nlm.nih.gov/pmc/articles/PMC3149148/?tool=pmcentrez

Discovering and Developing Molecules with Optimal Drug-Like Properties Editors: Allen C Templeton, Stephen R. Byrn, Roy J Haskell, Thomas E. Prisinzano

ISBN: 978-1-4939-1398-5 (Print) 978-1-4939-1399-2 (Online) 27 Sep 2014

http://link.springer.com/book/10.1007/978-1-4939-1399-2

> For More Trends, see

Dosing-Time Makes the Poison: Circadian Regulation and Pharmacotherapy Dallmann, Robert et al.

Trends in Molecular Medicine , Volume 22 , Issue 5 , 430-445 http://dx.doi.org/10.1016/j.molmed.2016.03.004

The right dose for every patient: a key step for precision medicine

Richard W. Peck

Nature Reviews Drug Discovery 15, 145-146 (2016) doi:10.1038/nrd.2015.22

http://www.nature.com/nrd/journal/v15/n3/full/nrd.2015.22.html

Animal-based studies will be essential for precision medicine

K. C. Kent Lloyd, Peter N. Robinson and Calum A. MacRae

Science Translational Medicine 17 Aug 2016:

Vol. 8, Issue 352, pp. 352ed12, DOI: 10.1126/scitranslmed.aaf5474

https://www.ncbi.nlm.nih.gov/pubmed/27535618

Can formulation and drug delivery reduce attrition during drug discovery and development—review of feasibility, benefits and challenges

S Basavaraj, Guru V. Betageri,

Acta Pharmaceutica Sinica B

Volume 4, Issue 1, February 2014, Pages 3-17

SI: Drug Delivery System and Pharmaceutical Technology

http://www.sciencedirect.com/science/article/pii/S2211383513001081

Circadian Timing in Cancer Treatments

Annual Review of Pharmacology and Toxicology Vol. 50: 377-421 (Volume publication date February 2010) DOI: 10.1146/annurev.pharmtox.48.113006.094626

http://www.annualreviews.org/doi/abs/10.1146/annurev.pharmtox.48.113006.094626?journalCode=pharmtox

Selected Further Reading II

(Miscellaneous Application Research Areas)

For comprehensive list, go to www.iprecio.com or contract Primetech Corporation.

> Drug Discovery: iPRECIO Micro infusion pumps referenced

Discovering and Developing Molecules with Optimal Drug-Like Properties

AAPS Advances in the Pharmaceutical Sciences Series, 2015

Editors: Templeton, A.C., Byrn, S.R., Haskell, R.J., Prisinzano, Th.E. (Eds.)

http://www.springer.com/jp/book/9781493913985?wt_mc=ThirdParty.SpringerLink.3.EPR653.About_eBook

Enabling Discovery Through Leveraging and Miniaturizing Pharmaceutical Principles and Processes

Discovering and Developing Molecules with Optimal Drug-Like Properties

AAPS Advances in the Pharmaceutical Sciences Series Volume 15, 2015, pp 95-140 Chapter 3, 27 Sep 2014,

http://link.springer.com/chapter/10.1007/978-1-4939-1399-2_3

Discovery Formulations: Approaches and Practices in Early Preclinical Development

Discovering and Developing Molecules with Optimal Drug-Like Properties

AAPS Advances in the Pharmaceutical Sciences Series Volume 15, 2015, pp 49-94 Chapter 2, 27 Sep 2014

http://link.springer.com/chapter/10.1007/978-1-4939-1399-2_2

Optimising in vivo pharmacology studies-Practical PKPD considerations

Journal of Pharmacological and Toxicological Methods Volume 61, Issue 2,

March-April 2010, Pages 146-156, Troubleshooting methods in pharmacology and toxicology

http://www.sciencedirect.com/science/article/pii/S1056871910000183

> Drug Discovery PK/PD

Translational pharmacokinetic-pharmacodynamic analysis in the pharmaceutical industry: an IQ Consortium PK-PD Discussion Group perspective

 $Drug\ Discov\ Today.\ 2017\ May\ 2.\ pii:\ S1359-6446(17)30002-8.\ doi:\ 10.1016/j.drudis.2017.04.015.$

https://doi.org/10.1016/j.drudis.2017.04.015

Implementation of pharmacokinetic and pharmacodynamic strategies in early research phases of drug discovery and development at Novartis Institute of Biomedical Research

Front. Pharmacol., 28 July 2014 | http://dx.doi.org/10.3389/fphar.2014.00174

http://journal.frontiersin.org/article/10.3389/fphar.2014.00174/abstract


Pharmacokinetics in Drug Discovery: An Exposure-Centred Approach to Optimising and Predicting Drug Efficacy and Safety.

Handb Exp Pharmacol. 2016;232:235-60. doi: 10.1007/164_2015_26.


http://link.springer.com/chapter/10.1007%2F164_2015_26

Example Pump implantation site and drug administration site


Intravenous Administration




SMP-310R



Intracerebral Administration



SMP-310R



Support Materials

Technical Note/Surgical Protocol:

- -Recommendation for Intravenous Administration.
- -Recommendations for Subcutaneous Administration.
- -Recommendations for Intraperitoneal Administration.
- -Recommendations for Intracerebral Administration.
- -Recommendations for Intrathecal Administration.

References

An Improved Method of Implanting a Programmable Continuous Infusion Pump in

68th AALAS National Meeting, October 15 to 19th 2017, Austin Convention Center, 500 E. Cesar Chavez Street, Austin, TX 78701, U.S.

Surgical Videos

Mouse Surgeries (SMP-300 / SMP-310R)



SMP-300 / SMP-310R with SC administration and general preparation video

https://drive.google.com/folderview?id=0B0pySJ1uXUqSVFBSVVAzTIZHaWc&usp=sharing



SMP-300 / SMP-310R with IP administration

https://drive.google.com/folderview?id=0B0pySJ1uXUqSd1BNdDVZeEFQUWM&usp=sharing



SMP-300 / SMP-310R with IV Jugular administration

https://drive.google.com/folderview?id=0B0pySJ1uXUqSbENyQ21nY2REcHM&usp=sharing



SMP-300 / SMP-310R with IV femoral administration

https://drive.google.com/folderview?id=0B0pySJ1uXUqSdlFtVFdJMFFncWM&usp=sharing



SMP-300 / SMP-310R with ICV administration

https://drive.google.com/folderview?id=0B0pySJ1uXUgSUGxHWkdLTXMwSEk&usp=sharing



Refilling Video and Refiling FAQ

https://drive.google.com/folderview?id=0B0pySJ1uXUqSX203d1I4bGsxOG8&usp=sharing

Rat Surgeries (SMP-200)



Surgery Training Videos

https://drive.google.com/folderview?id=0B0pySJ1uXUqSR2kzLVIMbWtRNUE&usp=sharing We have been working on surgical videos which we hope will help our users.

These are for the SMP-200 pumps you have been using

- > We have been working surgical videos unfortunately, they are not complete yet.
- > Feedback on the videos were provided by other surgeons (word document attached)
- > We have been working with Vetbiotech, www.vetbiotech.com to complete them.

Surgical Videos

From both Distributor and direct with Manufacturer (any way you want)

- · Phone
- · E-mail or fax
- Web meeting and training

iPRECIO® Key Features

- > Accurate patented Rotary Finger Method
 - Every pump is factory tested and calibrated
 - Better than ±5% accuracy
 - Programmable infusions protocols (simple and complex)
- > Totally implanted in subcutaneous space
- > Refillable (reservoir) percutaneously via refill port with re-sealable septum
- > With iPRECIO® catheters, test your drug's effects nearly anywhere
- > Easy to use software for infusion protocol programming



Implantable

The pump can be completely implanted in small laboratory animals subcutaneously. Thus, the animal moves freely without any restrain (i.e. tethering) during drug infusion. Additionally, infection risk is reduced, and the animal is likely to be significantly less stressed than in a tethered infusion model.



Refillable

You can replenish or exchange saline and/or any medical fluid in the pump via percutaneous access to the pump refill septum and reservoir after implantation of the pump. Recovery from surgery or washouts may be planned with saline in the reservoir. Long-term drug infusion can be maximized to battery life of the pump.




Precision


The technology driving the infusion is a patented "Rotary Finger" method. This method is a unique form of peristalsis. The precise "micro-stick" pushes a rubber tube in the pump in a uniform and sequential manner. The accuracy of iPRECIO is +/-5%.

Programmable


> SMP-310R 15 steps for flow rate or dose programming: 0.0-10.0 ul/hr with repeat mode


Each flow profile may contain up to 15 doses or flow rate steps. A single step would mean a fixed continuous dose or flow-rate for the study duration. A more complex infusion profile will contain more than 1 step and may contain up to 15 steps. KVO and dead volume flushing functions may be programmed within the 15 programmable steps.







> SMP-200 10 steps for flow rate programming: 0, 0.2, 0.5 & 1.0 - 30.0 ul/hr with repeat mode










iPRECIO® is an Ultimate Choice

This implantable infusion pump uses a patented, microprocessor controlled peristalsis mechanism for accurate controlled flow. It is the only implantable and programmable pump for small laboratory animals. iPRECIO® can infuse fluids continuously for as long as six months and it can be refilled via a percutaneously accessible port.




iPRECIO® Pump's Structure


> SMP-310R





> SMP-200





iPRECIO® Management System

> SMP-310R





iPRECIO® Management System is sold as IMS-310R which consists of data communication device (UCD-X10R) and Management Software, User Manual.

> SMP-200



iPRECIO® Management System consists of:

- Data Communication Device
- USB cable, 2 AAA batteries
- iPRECIO® Management Software Installation CD
- iPRECIO® User Manual

iPRECIO® Battery Life

> SMP-310R

Com.	Per n	ninute	Every 2	2 hours	Every 6	6 hours	Every 2	4 hours	No	ne
Flow	Driving									
rate	hours	days								
0.1	157	6.5	528	22.0	1063	44.3	1542	64.3	1628	67.8
0.5	155	6.5	476	19.8	887	37.0	1214	50.6	1266	52.8
1.0	153	6.4	428	17.8	742	30.9	959	40.0	991	41.3
5.0	137	5.7	263	11.0	344	14.3	357	14.9	362	15.1
8.0	127	5.3	207	8.6	243	10.1	243	10.1	245	10.2
10.0	121	5.0	178	7.4	196	8.2	200	8.3	201	8.4

Flow Rate Unit: µL/hr

> SMP-200

Flam Data	Infu	Tatal Maluma	
Flow Rate	Time (h)	Days (approx.)	Total Volume
30.0 μL/hr	196 hr	1 week	5.8 ml
19.0 μL/hr	307 hr	1.8 weeks	5.8 ml
8.5 µL/hr	669 hr	1 month	5.6 ml
1.0 µL/hr	4,328 hr	6 months	4.3 ml

^{*} Table above outlines the maximum battery life for the programmed protocol and pump switch on time. Exact battery life will be dependent on pump switch on time, programmed infusion protocol, and selected communication availability(Com.). iPRECIO Management software helps the user calculate battery life for selected programming.

Model	SMP-310R IMS-310R (new!!)	SMP-300 IMS-300 (discontinued)	SMP-200 IMS-200	
Appearance of the pump	24.8(L) x 15.0(W) x 7.2 (H) mm	24.8(L) x 15.0(W) x 7.2 (H) mm	Control of the contro	
	Max. height 7.5mm	Max. height 7.4mm	38.7 (L) X 19.2 (W) X 9.7 (H) mm	
Туре	Implantable SC	Implantable SC	Implantable SC	
Volume / Weight	2.26cc / 3.4g	2.15cc / 3.3g	7.20cc / 7.9g	
Animal Species	Mouse or larger	Mouse or larger	Rats or larger	
Reservoir Volume	130 µL	130 µL	900 μL	
Flow Rate (Setting Resolution)	0.0 – 10.0 μL/hr (0.1μL/hr)	0.0 – 10.0 μL/hr (0.1μL/hr)	0.0, 0.2, 0.5&1.0 – 30.0μL/hr (0.1μL/hr)	
Flow Steps / Repeat	15 / Yes	15 / Yes	10 / Yes	
	0 & 0.1 ul/hour 67 days	0 & 0.1 ul/hour 46 days	0, 0.2, 0.5, 1 µl/hour - 6 mths	
Battery Life	1 μl/hour up to 41 days	1 μl/hour up to 33 days	2.5 µl/hour - 86 days	
	10 μl/hour up to 8 days	10 μl/hour up to 9 days	30 μl/hour - 8 days	
Programmable	Wireless Preprogrammable	Wireless Preprogrammable	Preprogrammed prior to implantation	
Wireless Distance	1 - 6m	< 1m	-	
Communication Availability	1m, 1h, 2h, 4h, 6h, 12h, 24h and NONE (8 choices)	1m, 2h, 6h, 24h and None (4 choices)	-	

Compatible solvents for SMP-300, 310R and SMP-200

- * Tested for both SMP-200 & SMP-300 / SMP-31R
- * Tested in SMP-200 Pump Only

(same materials and manufacturing process) and expected to be compatible when compatible. Also, not compatible when not compatible.

Compatible Solvents

Acids, with pH 2 or weaker *

Bases, with pH less than 13 *

Buffered Phosphate Saline (PBS) *

Culture Media (1% benzyl alcohol) *

Cyclodextrin *

Dextrose, up to 5% in water or saline *

N,N-Dimethyl formamide (DMF), up to 25% in water *

DMSO 50% and water or saline 50% *

DMSO, up to 50% in ethanol (≤15%) and water *

DMSO 5% and PEG400 95%

50% DMSO + 50% Propylene Glycol *

DMSO 50% and water 50% *

DMSO 50% + 15% ethanol and 35% water *

Dulbecco's Modified Eagle Medium (D-MEM) (1X), liquid *

Ethanol, up to 50% in water *

Glycerin, up to 75% in water *

Glycerol 100% *
1-Methyl-2-Pyrrolidone, up to 12.5% in water *

Propylene Glycol *

Ringer's solution (without lactate) *

Saline, 0.9% (or other aqueous salt solution) *

Triacetin, up to 5% in water *

Tween 80, up to 2% in water *

Water, distilled *

PEG200 100% *

Solutol® 15% in water *

Viscosity up to 20 cp is ok.

(Higher viscosity not tested due to the use of 27G needles.

Difficulty to aspirate solution with 27G needle)

Short term use only (1 - 2month)

PEG300 100% * (< 45 days)

PEG400 100% *

Cremophor EL 25% in water * (< 30 days)

PEG400/Propylene Glycol/Water 30 : 50 : 20 * (< 30 days)



www.iprecio.com E-mail : iprecio@primetech.co.jp

a division of Harvard Bioscience, Inc.

Rev16 Apr. 2020 Distributed and Supported by: