

Instructions to update Ponemah for compatibility with HD-S20 and HD-S1

To use the new HD-S20 and/or HD-S1 implants in Ponemah versions v6.33 and v5.20-SP9, an updated product definition file is required for the DSI software. Ponemah must be at least Ponemah version 6.30 or 5.20-SP8 to install this new file.

This updated product definition file can be obtained from the DSI website and can be upgraded on your PC(s) using the following steps:

1. Open a web-browser and navigate to the following link:
<https://www.datasci.com/resources/software-and-firmware-updates>
2. Download HD-S20 and HD-S1 Product Definitions.zip and extract the files
3. Close Ponemah if open
4. Open Windows File Explorer on your PC and navigate to the following folder:

C:\ProgramData\DSI\DACSS

5. Find the file **EDemProductDefinition.xml** and rename the file to **EDemProductDefinition.xml.OLD**
6. Copy the **EDemProductDefinition.xml** to the **C:\ProgramData\DSI\DACSS** folder

The next time you start your DSI software, you will be able to successfully find, configure, and acquire with your HD-S20 and/or HD-S1 implants.

If you are still unable to use your new implants please contact DSI Technical Support at:
Support@datasci.com