



PhysioDetect

Work Smarter, Not Harder.

PhysioDetect helps you confidently expedite your study setup, and ongoing study management, with immediate feedback regarding critical PhysioTelTM Digital large animal implant information.

User-friendly Dashboard

- View implants in range with the ability to filter only implants or frequencies of interest
- Access critical implant information
- Targeted search to find implants by serial number

Simplified Device Detection

- Confirm implant activation
- View implant frequency
- Save implants to an active "watch list"

Flexible Subscription Service

- Instant version updates
- No setup required
- Affordable annual contracts
- Preconfigured hardware to ensure security and compliance



Confirm serial number, implant mode, battery life, and frequency setting from the PhysioDetect app.

DEVICE	STATUS	BATTERY	
00005921 Test Subject 10-1	ACTIVE	75%	A1-F1
00005922 Test Subject 10-2	ACTIVE	98%	A1-F1
00005923 Test Subject 10-3	WARNING	5%	A1-F1
00005924 Test Subject 10-4	ACTIVE	22%	A1-F1
00005925 Test Subject 10-5	STANDBY	?	A1-F1
00005926 Test Subject 10-6	ACTIVE	55%	A1-F1
00005931 Test Subject 11-1	SCANNING	?	B1-F1
00005932 Test Subject 11-2	SCANNING	?	B1-F1
00005933 Test Subject 11-3	SCANNING	?	B1-F1
00005934 Test Subject 11-4	SCANNING	?	B1-F1

Contact your DSI representative today to learn more about PhysioDetect.