

## Recommended Instruments For Large Animal Surgery

The surgical instruments, along with their part numbers from Fine Science Tools, are listed below.

Part Number	Description
11006-12	Adson Forceps-straight, serrated
11027-12	Adson-Brown Tissue Forceps-straight, with teeth
13019-14	Kelly Hemostat Forceps, curved
13018-14	Kelly Hemostat Forceps, straight
13009-12	Halsted Mosquito Forceps, curved
13008-12	Halsted Mosquito Forceps, straight
11617-12	Debaquey Forceps
14010-17	Mayo Scissors-straight, 15 cm
10019-14	Metzenbaum Scissors-curved, 14.5 cm
12002-14	Olsen-Hegar Needle Holder
11095-09	Backhaus Towel Clamps

A hollow trocar or skin tunneling needle is also necessary to tunnel the pressure catheter or ECG leads subcutaneously from the transmitter site to the intended placements. They are available from Alice King Chatham Medical Arts (AKCMA) or Chiron Bioscience.

Fine Science Tools, Inc.

Telephone: (1-800) 521-2109 or (1-650) 349-1636

Fax: (1-800) 523-2109 or (1-650) 349-3729

Website: [www.finescience.com](http://www.finescience.com) (A list of offices in other countries can also be found here.)

Alice King Chatham Medical Arts

Telephone: (1-310) 970-1063

Fax: (1-310) 970-0121

Website: <http://www.akcma.com>

Chiron Bioscience Limited

Website: [www.chironbioscience.com](http://www.chironbioscience.com)

For a description and ordering information, search Chiron Skin Trocar and cannula

A wound clip applier with wound clips can be purchased from Fisher Scientific (Part Number NC9154268).

Fisher Scientific

Telephone: (1-800) 766-7000

Fax: (1-800) 926-1166

Website: [www.fishersci.com](http://www.fishersci.com)